

TC WIRE CUTTING PLIER; ANGLED

soft wire up to Ø2.0mm; hard to Ø1.7mm

double convertible; 17.0cm (6.8")


SKU: 190-773-17

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX159
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-523-16
Medicon	14.30.17



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX159
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-523-16

Medicon

14.30.17