

## TC WIRE CUTTING PLIER; CUTTING WIDTH 18M

soft wire up to Ø2.0mm; hard to Ø1.2mm

double convertible; 15.0cm (6")



**SKU:** 190-770-15

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Forceps</a>
<b>Aesculap</b>	LX150
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-524-15
<b>Medicon</b>	14.30.65

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire Forceps</a>
<b>Aesculap</b>	LX150
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-524-15

**Medicon**

14.30.65