

TC WIRE CUTTING PLIER

soft wire up to Ø3.0mm; hard to Ø2.2mm

double convertible; 23.0cm (9.2")


SKU: 190-774-23

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX156
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-523-24
Medicon	76.45.73



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Wire Forceps
Aesculap	LX156
MATERIAL	Tungsten Carbide(TC)
KLS Martin	22-523-24

Medicon

76.45.73