

TC WIRE CUTTING PLIER; CUTTING WIDTH 10MM

soft wire up to Ø2.0mm; hard to Ø1.5mm

double convertible; 18.0cm (7.2")

SKU: 190-772-18

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Wire Forceps |
| Aesculap | LX158 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 22-526-18 |
| Medicon | 14.30.28 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Wire Forceps |
| Aesculap | LX158 |
| MATERIAL | Tungsten Carbide(TC) |
| KLS Martin | 22-526-18 |

Medicon

14.30.28