

## TC WIRE TWISTER

straight; serrated

18.5cm (7.4")




**SKU:** 051-242-18

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire twisting Forceps</a>
<b>Aesculap</b>	BM113
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>
<b>KLS Martin</b>	22-506-18
<b>Medicon</b>	14.44.67
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Wire twisting Forceps</a>
<b>Aesculap</b>	BM113
<b>MATERIAL</b>	<a href="#">Tungsten Carbide(TC)</a>

**KLS Martin**

22-506-18

**Medicon**

14.44.67

**SPECIFICATIONS**

[serrated](#), [straight](#)