

TELESCOPE HD+2 FOR LAPAROSCOPY

HYSTEROSCOPY, CYSTOSCOPY

Ø 2.9mm; 12°; 300 mm; autoclavable



SKU: 601-103-12HD2

Categories: [Endoscopy](#), [Equipment](#), [Laparoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)



12°

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [GUBBINI/Gyn](#), [Gynecology](#), [Laparoscopy](#), [Urology](#)