

TENDON INTERLACING FORCEPS

straight

15.0cm (6")



SKU: 195-234-15
Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	FO950
KLS Martin	23-682-15
Medicon	78.12.16
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Seizing Instruments
Aesculap	FO950
KLS Martin	23-682-15
Medicon	78.12.16

SPECIFICATIONS

[straight](#)