

TENDON STRIPPER

Ø 4.0mm; flexible

23.0cm (9.2")

SKU: 195-265-04

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
KLS Martin	23-654-40
Medicon	78.20.54

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
KLS Martin	23-654-40
Medicon	78.20.54