

TERRIER HAEMOSTATIC FORCEPS

straight; serrated; 1X2 teeth

8.5cm (3.4")



SKU: 040-016-08

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
KLS Martin	13-154-08
AUTHORS	TERRIER
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
KLS Martin	13-154-08
AUTHORS	TERRIER
SPECIFICATIONS	serrated , straight