

TOENNIS NEEDLE HOLDER

straight; serrated; without catch

18.0cm (7.2")

SKU: 050-362-18

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM510
KLS Martin	20-574-18
AUTHORS	TOENNIS
Medicon	11.34.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM510
KLS Martin	20-574-18

AUTHORS

[TOENNIS](#)

Medicon

11.34.18

SPECIFICATIONS

[serrated](#), [straight](#)