

TONG TYPE ELECTRODE BIPOLAR; 3MM; INSERT

for bipolar coagulation forceps

628-101-00 and 628-105-00; 330mm

No CE mark! No labeling!

SKU: 628-140-53

Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes