

## TOURNIQUET LIGATURE CATCHER

Ø2.5mm; 26.0cm




**SKU:** 141-028-25

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Ligature Catchers, Tourniquets](#)  
**KLS Martin** 24-550-25



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Ligature Catchers, Tourniquets](#)  
**KLS Martin** 24-550-25