

## TOURNIQUET LIGATURE CATCHER

Ø2.0mm; 26.0cm


**SKU:** 141-028-20

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Catchers, Tourniquets</a>
<b>KLS Martin</b>	24-550-20



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Catchers, Tourniquets</a>
<b>KLS Martin</b>	24-550-20