

**TOWEL FORCEPS FOR PAPER CLOTHS**


curved; ball/socket

13.0cm (5.2")



**SKU:** 030-024-13

**Categories:** [Open Surgery](#), [Hemostats](#)



**ZUSÄTZLICHE INFORMATIONEN**

|                       |                               |
|-----------------------|-------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>          |
| <b>TYPE</b>           | <a href="#">Towel Forceps</a> |
| <b>Aesculap</b>       | BF464                         |
| <b>KLS Martin</b>     | 13-924-14                     |
| <b>Medicon</b>        | 17.55.23                      |
| <b>SPECIFICATIONS</b> | <a href="#">curved</a>        |

**ZUSÄTZLICHE INFORMATIONEN**

|                   |                               |
|-------------------|-------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>          |
| <b>TYPE</b>       | <a href="#">Towel Forceps</a> |
| <b>Aesculap</b>   | BF464                         |
| <b>KLS Martin</b> | 13-924-14                     |
| <b>Medicon</b>    | 17.55.23                      |

## SPECIFICATIONS

[curved](#)