

TOWEL FORCEPS FOR PAPER CLOTHS

curved; ball/socket

11.5cm (4.5")




SKU: 030-024-11

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF465
KLS Martin	13-924-11
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
Aesculap	BF465
KLS Martin	13-924-11
SPECIFICATIONS	curved