

## TOWNLEY CALIPER

11.5cm (4,6")



**SKU:** 100-266-10

**Categories:** [Open Surgery](#), [Diagnostics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Measuring instruments</a>
<b>Aesculap</b>	AA847
<b>KLS Martin</b>	17-409-01
<b>AUTHORS</b>	<a href="#">TOWNLEY</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Measuring instruments</a>
<b>Aesculap</b>	AA847
<b>KLS Martin</b>	17-409-01
<b>AUTHORS</b>	<a href="#">TOWNLEY</a>