

## TRIPOD POLYP RETRIEVER WITHOUT HANDLE;

teflon tube, with 3 sharp hooks

1.8mm, length: 180cm



**SKU:** PK-7GZPK1840180C

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Foreign Body Forceps](#)