

TRIPOD POLYP RETRIEVER

teflon tube, with 3 round

2.3mm, length 250cm; single use



SKU: PK-7GREW2340250C-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopia](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopia](#)