

TROCAR FOR HASSON, CONICAL BLUNT

Ø 11.0mm; NL:150mm

for sleeve 626-121-99

SKU: 626-157-99
Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Trocar Systems
TYPE	Trocars

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Trocar Systems
TYPE	Trocars