

TUBE ROLLER SQUEEZING FORCEPS

teflon roller, aluminium blue anodised

16.0cm (6.4")



SKU: 032-210-00
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
SPECIFICATIONS	straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
SPECIFICATIONS	straight