

TUBING CLAMP

straight; serrated

18.0cm (7.2")




SKU: 032-002-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN912R
KLS Martin	13-944-18
Medicon	17.72.78
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN912R
KLS Martin	13-944-18
Medicon	17.72.78

SPECIFICATIONS

[serrated](#), [straight](#)