

TUBING CLAMP

straight; serrated

18.0cm (7.2")



SKU: 032-002-18

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Tubing Clamps](#)

Aesculap AN912R

KLS Martin 13-944-18

Medicon 17.72.78

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

Aesculap AN912R

KLS Martin 13-944-18

Medicon 17.72.78

SPECIFICATIONS

[serrated](#), [straight](#)