

TUBING CLAMP

straight; serrated

20.0cm (8")



SKU: 032-002-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Tubing Clamps](#)

Aesculap AN913

KLS Martin 13-944-20

Medicon 17.72.80

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

Aesculap AN913

KLS Martin 13-944-20

Medicon 17.72.80

SPECIFICATIONS

[serrated](#), [straight](#)