

TUBING CLAMP

straight; serrated; with Clip

18.0cm (7.2")



SKU: 032-004-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Tubing Clamps](#)

Aesculap AN912

KLS Martin 13-946-18

Medicon 17.72.79

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

Aesculap AN912

KLS Martin 13-946-18

Medicon 17.72.79

SPECIFICATIONS

[serrated](#), [straight](#)