

## TUBING CLAMP

straight; serrated; with Clip

18.0cm (7.2")



**SKU:** 032-004-18

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Tubing Clamps](#)

**Aesculap** AN912

**KLS Martin** 13-946-18

**Medicon** 17.72.79

**SPECIFICATIONS** [serrated](#), [straight](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Tubing Clamps](#)

**Aesculap** AN912

**KLS Martin** 13-946-18

**Medicon** 17.72.79

## SPECIFICATIONS

[serrated](#), [straight](#)