

TUBING CLAMP

straight; serrated

15.5cm (6.2")



SKU: 032-002-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

Aesculap AN910

KLS Martin 13-944-15

Medicon 17.72.75

SPECIFICATIONS [serrated](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tubing Clamps](#)

Aesculap AN910

KLS Martin 13-944-15

Medicon 17.72.75

SPECIFICATIONS

[serrated](#), [straight](#)