

TWIST DRILL FOR HAND DRILL

Ø 2.5mm

SKU: 190-261-25
Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR041
KLS Martin	22-145-01
Medicon	76.50.11

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR041
KLS Martin	22-145-01
Medicon	76.50.11