

TWIST DRILL FOR HAND DRILL

Ø 4.0mm

SKU: 190-261-40
Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR044
KLS Martin	22-145-04
Medicon	76.50.14

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
Aesculap	FR044
KLS Martin	22-145-04
Medicon	76.50.14