

## TYDING TONSIL SNARE

with one straight tip



**SKU:** 132-524-01

**Categories:** [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tonsil Snares</a>
<b>Aesculap</b>	OM781
<b>KLS Martin</b>	35-525-01
<b>AUTHORS</b>	<a href="#">TYDING</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tonsil Snares</a>
<b>Aesculap</b>	OM781
<b>KLS Martin</b>	35-525-01
<b>AUTHORS</b>	<a href="#">TYDING</a>