

VAPORIZATION ELECTRODE

Ø 3mm; toothed; 24 / 26Fr.
 yellow/blue; bipolar system



SKU: 630-240-45B
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes