

VAPORIZATION ELECTRODE

Ø 5mm; toothed; 24 / 26Fr.

yellow; monopolar



SKU: 630-240-50
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes