

## VERES INSUFFLATION CANNULA

Ø 2.0mm; FL 150mm

stopcock; Luer-Lock; stainless steel



**SKU:** 618-550-15

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [VERRES Cannulas](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [VERRES Cannulas](#)