

## VERES INSUFFLATION CANNULA; HIGH FLOW

Ø 2.7mm; FL 120mm

stopcock; Luer-Lock; stainless steel



**SKU:** 618-555-12

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [VERRES Cannulas](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [VERRES Cannulas](#)