

WAGENER EAR HOOK

2.0mm

14.0cm (5.6")

SKU: 111-588-05

Categories: [Open Surgery](#), [Otology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Loops, Hooklets
KLS Martin	36-585-05
AUTHORS	WAGENER
Medicon	62.16.01



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Loops, Hooklets
KLS Martin	36-585-05
AUTHORS	WAGENER
Medicon	62.16.01