

WANGENSTEEN NEEDLE HOLDER

straight; serrated

27.0cm (10.8")



SKU: 050-160-27

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM272
KLS Martin	20-420-27
AUTHORS	WANGENSTEEN
Medicon	10.40.28
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM272
KLS Martin	20-420-27

AUTHORS

[WANGENSTEEN](#)

Medicon

10.40.28

SPECIFICATIONS

[serrated](#), [straight](#)