

## WASHING TUBE

standard tip

2.3mm, length 250cm; single use



**SKU:** PK-8MSEW2312250-S

**Category:** [Endoscopy](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)