

WASSMUND RETRACTOR

2 teeth; blunt; 33X21mm

20.5cm (8.2")


SKU: 064-027-02

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	DO409
KLS Martin	15-297-21
AUTHORS	WASSMUND
Medicon	68.18.91
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Retractors
Aesculap	DO409
KLS Martin	15-297-21

AUTHORS

[WASSMUND](#)

Medicon

68.18.91

SPECIFICATIONS

[blunt](#), [toothed](#)