

WAUGH DRESSING FORCEPS

straight; serrated

15.0cm (6"); width of jaw 1.8 +/-0.2

SKU: 020-130-15
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD154
AUTHORS	WAUGH
Medicon	06.42.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD154
AUTHORS	WAUGH
Medicon	06.42.15

SPECIFICATIONS

[serrated](#), [straight](#)