


WEITLANER SELFRETAINING RETRACTOR

sharp; 3X4 teeth
 20.0cm (8")



SKU: 070-670-20
Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV200
KLS Martin	15-734-20
AUTHORS	WEITLANER
Medicon	24.22.20
SPECIFICATIONS	sharp , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV200
KLS Martin	15-734-20

AUTHORS

[WEITLANER](#)

Medicon

24.22.20

SPECIFICATIONS

[sharp](#), [toothed](#)