


## WILDCAT GUIDE WIRE,

Ø 0.025", length: 480cm (single use)

teflon coated; platin tip J style



**SKU:** PK-3LQPK0550480XJ-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)