

WIRE LOOP FOR TONSIL SNARE

Ø 50 x 0.4mm; package of 100 pieces




SKU: 132-520-50

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Snares
Aesculap	OM758
KLS Martin	39-511-50



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Snares
Aesculap	OM758
KLS Martin	39-511-50