

## WORKING ELEMENT

active cutting; closed handle

for Ø 4mm telescopes;HF-connection above

Only for monopolar use



**SKU:** 630-235-11

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Resectoscopy 4.0](#)  
**TYPE** [Working Elements](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**RANGE** [Resectoscopy 4.0](#)  
**TYPE** [Working Elements](#)