

WORKING ELEMENT

passive cutting; closed handle
 for Ø 4mm telescopes; HF-connection down
 Only for monopolar use



SKU: 630-235-60
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements