

WORKING ELEMENT

passive cutting; closed handle
 for Ø 4mm telescopes;HF-connection above
 Only for monopolar use



SKU: 630-235-61
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements