


**WORKING ELEMENT; BIPOLAR; TITAN**

active cutting; closed handle  
 WO compatible / GYRUS PK  
 incl. 1 cone sleeve / Titanium blue



**SKU:** 630-351-05  
**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

**ZUSÄTZLICHE INFORMATIONEN**

|       |                                  |
|-------|----------------------------------|
| UNIT  | <a href="#">pcs.</a>             |
| RANGE | <a href="#">Resectoscopy 4.0</a> |
| TYPE  | <a href="#">Working Elements</a> |

**ZUSÄTZLICHE INFORMATIONEN**

|       |                                  |
|-------|----------------------------------|
| UNIT  | <a href="#">pcs.</a>             |
| RANGE | <a href="#">Resectoscopy 4.0</a> |
| TYPE  | <a href="#">Working Elements</a> |