


WORKING ELEMENT PASSIVE TITANIUM BLUE

KS compatible incl. 2 Cone sleeves
for 4mm/302mm Telescope

GYRUS PK



SKU: 630-350-55
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements