

## YANKAUER SUCTION TUBE

removable suction tip; Luer connector

inner Ø 2.0mm; outer Ø 6.0mm; 29.0cm

**SKU:** 085-000-27

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF944
<b>KLS Martin</b>	18-543-27
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	09.85.01

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF944
<b>KLS Martin</b>	18-543-27
<b>AUTHORS</b>	<a href="#">YASARGIL</a>

**Medicon**

09.85.01