

YASARGIL APPLYING FORCEPS

bayonet; for Standard vessel clips

17.0cm (6.8")



SKU: 244-601-17
 Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
Aesculap	FE580
AUTHORS	YASARGIL
Medicon	58.54.97
SPECIFICATIONS	jaw rigid

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
Aesculap	FE580
AUTHORS	YASARGIL
Medicon	58.54.97

SPECIFICATIONS

[jaw rigid](#)