

YASARGIL APPLYING FORCEPS; JAW ROTATEABLE

bayonet; for standard vessel clips

21.0cm (8 1/4")

SKU: 244-603-21
Categories: [Open Surgery](#), [Neurosurgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotateable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotateable