

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 13.0cm (5.2")rotatable
 for Mini vessel clips



SKU: 244-680-13
Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
Medicon	58.37.94
SPECIFICATIONS	jaw rigid , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
Medicon	58.37.94
SPECIFICATIONS	jaw rigid , shaft rotatable